

# sigInspect User's Manual

Eduard Bakstein, [eda@zzz.cz](mailto:eda@zzz.cz)

December 12, 2016

**sigInspect** is a Matlab GUI tool for inspecting multi-channel recordings. Is particularly useful for extracellular microelectrode recordings and allows *visualization*, *playback* and *annotation* of segments containing *artifacts*. SigInspect was developed by the NEURO group (<http://neuro.felk.cvut.cz>) at the Department of Cybernetics, Faculty of Electrical Engineering, Czech Technical University in Prague. SigInspect repository with the most recent version can be found at <https://github.com/ebakstein/sigInspect>.

## Contents

<b>Overview</b>	<b>1</b>
<b>Quick start</b>	<b>2</b>
Signal . . . . .	2
Loading a signal - the simple way . . . . .	2
Loading signals from a mat-file . . . . .	2
Using the GUI . . . . .	2
Settings . . . . .	4
<b>Advanced data loading - sigInspectDataInterface</b>	<b>4</b>
In OOP terms . . . . .	6
Reference . . . . .	6
Example interface . . . . .	7
Default interface: sigInspectDataBasic . . . . .	7
<b>Automatic artifact annotation (Beta)</b>	<b>8</b>
Artifact types handling . . . . .	8
Artifact classification methods . . . . .	8
<b>About + license</b>	<b>9</b>
<b>FAQ</b>	<b>9</b>

## Overview

**sigInspect** is a graphical user interface (GUI) application for Matlab, developed for inspection and annotation of extracellular microelectrode recordings (MER). The tool allows concurrent visualization of multiple parallel channels, playback and spectrogram plot of selected channels (to identify specific firing patterns) and especially annotation of individual channels. Parallel signals are displayed in one-second segments, which can be easily annotated using the GUI controls or shortcut keys. SigInspect also contains tools for **automatic annotation of microelectrode recording data**, based on published methods. This can provide starting point for manual annotation and speed up the data-cleaning process considerably.

## Quick start

### Signal

A signal in `sigInspect` is assumed to be a **row vector** with samples in columns, a **multi-channel signal is a matrix with channels in rows**. `SigInspect` shows always one second segment of all channels in current signal. A signal can have one to unlimited number of channels.

### Loading a signal - the simple way

The simplest way to start viewing and annotation is to call `sigInspect` with signals as a parameter. Multiple signals can be passed in a cell array instead - user can then choose the displayed signal using the "signal" selector. The second parameter is sampling frequency in Hz (default value:  $24kHz$  is set if second parameter is omitted)

```
1 % 1 - single multi-channel signal (matrix as input)
2 size(signal)                % C x N matrix (C = channels, N = samples )
3 sigInspect(signal, samplingFreq); % signal: chan. in rows, samples in columns
4
5 % 2 - multiple signals (cell array as input)
6 s={signal1,signal2,signal3};
7 sigInspect(s, samplingFreq);
```

### Loading signals from a mat-file

In case `sigInspect` is called without parameters, an *Open file dialog* will pop up upon initialization, asking you to select a `*.mat` file with data. The mat-file has to contain signals in a cell array - the format described above - in a **variable called signal, signals or data**. Otherwise, your signals will not be found and the application will issue an error.




Figure 1: Main window of the sigInspect GUI with loaded 4-channel signal

### Using the GUI

The sigInspect GUI with loaded signals can be seen in Figure 1. The main window consists of *time series plot* with selected second of given (multi-channel) signal and *spectrogram*, calculated from currently

selected channels.

## Basic controls

**Switching between seconds** can be done using the  $\leftarrow$  left and right  $\rightarrow$  arrow keys. When you reach the last second of current signal and hit right, you get to the first second of subsequent signal and vice versa. You can of course use corresponding buttons called "previous" and "next" or "signal" and "second" selector boxes with the same effect.

To have a clear view of all your channels, you will probably need to **adjust gain** (amplitude) of the signals. This can be done either by the "gain" slider on the left or by pressing  $\uparrow$  up or  $\downarrow$  down arrow keys.

## Channel selection

There are several operations that affect only the selected channels A) annotation (when you mark a specific 1s signal segment as an artifact) B) audio playback (when audio is on) and C) spectrogram calculation. **Channel selection** can be done using the numpad keys (1-9), number keys (top row of your keyboard with ENG layout)<sup>1</sup> or with checkboxes on the right, next to each signal. You can also **invert selection** using the "i" key or **select all** with the "a" key.

## Audio

SigInspect allows **playback of thresholded signals**, which is useful for artifact identification, as well as for simple overview of neuronal activity<sup>2</sup>. To hear audio playback, "Sound on" checkbox at the bottom has to be checked (use "o" shortcut to toggle sound), at least one channel must be selected and threshold has to be set lower than at least some parts of some of the selected channels. **Threshold** is indicated by the black horizontal lines surrounding each channel in the time series plot and can be adjusted by the "thr" slider on the left or by hitting the + and - keys. SigInspect **plays only signal parts exceeding currently set threshold in selected channels**. This can be used to play back only firing patterns of close neurons. If you wish to hear the whole signal including neuronal background, set threshold to 0;

If you wish to **replay** current second, use the "r" key, hit "Replay" button or change channel selection.

## Annotation

Annotation can be done by the artifact type buttons at the bottom or corresponding "F1"- "F12" keys. The number and types of artifacts are up to you and correspond to the ARTIFACT\_TYPES field in the settings. **When you hit a button, you assign it to all currently selected channels**. If the currently selected channels have this artifact type already assigned, you will unmark them. **You can assign from zero to all your defined artifact types to each second of each channel**. The labels on each button also show numbers of all selected channels, that have given artifact assigned.

Annotation can be **saved to a \*.mat file** either by clicking the save icon in the toolbar, the "Save all button at the bottom" or by hitting "s". By default, annotation is saved to current directory as sigInspectAnnotyyyy-mm-dd-HHMMSS.mat, where yyyy-mm-dd-HHMMSS is current date and time (to change this, see ANNOT\_DEFAULT\_FILENAME in the settings). Annotation can be loaded from a mat-file using the load icon in the toolbar.

The **format of saved annotation** is a cell array of the same length as the number of signals. Each cell contains annotation for given multichannel signal in a  $C \times S \times A$  logical matrix, where  $C$  is the number of channels,  $S$  signal length in seconds and  $A$  number of defined artifact types. For example, `annotation(3,2,1);` is one when the 3rd channel contains artifact of type 1 in 2nd second. Saved annotation also contains interface type, artifact types and signalIds to verify that you are loading annotation that fits to currently loaded data (when metadata do not fit, sigInspect offers you to match/repair the metadata).

If you switch the overview window on, you can also show signal segments annotated as artifact in color.

<sup>1</sup>Number shortcuts work obviously only up to 10 channels (tenth channel is controlled by 0)

<sup>2</sup>This method is commonly used in Deep Brain Stimulation surgery for nuclei identification

## Overview window

By checking the "overview checkbox" or by typing "v" you turn on the overview window, which **shows whole current multi-channel signal**. Current second is marked by red rectangle - you can select any second by clicking on it in the overview window.




Figure 2: Overview window showing 4-channel signal from Figure 1 with annotation

The overview window **can also display annotated artifacts** - just select the "all" radio button at the top. All clean segments with no artifact assigned will remain in default color (blue), while artifact seconds will turn red (can be changed in the settings).

To see how long is the current **longest artifact-free segment**, select maxCleanSegment radio button. This will mark all shorter clean segment red, so that you have a quick idea how long clean contiguous segment of your signal is left.

The overview window uses further resampling to speed up display (see OVERVIEW\_DECIMATE\_FACTOR in the Settings section of this guide)

## Keyboard shortcuts

For the sake of speed and usability, many controls can be accessed through simple keyboard shortcuts (mostly single-button). Most of the shortcuts are written directly in the GUI: you can e.g. press "o" to toggle sound on/off, "f" to toggle spectrogram or hit number keys to select or deselect channels.

## Settings

To accomodate sigInspect to your needs, it may be necessary to alter some of its default settings. If you use your own data interface, the settings can be changed easily from there (see the "Advanced data loading" Section for more details). If you use sigInspect without defining your own sigInspectDataInterface, you can change settings easily using example code in Section *Default interface*.

All the main available settings are summarised in Table 1 - the most important ones are in bold. Settings are being continuously added to the sigInspect, see the source code of sigInspect / ... sigInspect.OpeningFcn for a full list.

## Advanced data loading - sigInspectDataInterface

To provide the possibility to connect sigInspect to your existing database and to make configuration for multiple data sources easy, you can implement the sigInspectDataInterface and connect the GUI to whatever data source you are using. When you call sigInspect, you simply provide an instance of your implemented interface as the only parameter.

The two main things you have to implement are two methods: getSignalIds() which returns a list of signal identifiers (a cell array of strings) and getSignalsById(signalId), which returns a (multi)channel signal matrix for given signalId. You can also override the default settings within your interface - this way, you can easily adapt sigInspect for multiple data sources you work with. Even

Table 1: Overview of available sigInspect settings

field name	default value	meaning
basics & signal handling		
<b>SAMPLING_FREQ</b>	24000	signal sampling frequency in Hz <sup>3</sup>
<b>PLOT_STEP</b>	150	distance on the y axis between channels
<b>PLOT_CHANNELS</b>	5	max. number of prallel channels to plot
DECIMATE_FACTOR	1	decimate signals by this factor upon load (e.g. to save memory & processing time)
ADAPT_GAIN_TO_SIGNAL	1	automatically adapt gain slider to signal
ADAPT_GAIN_QUANTILE	.002	use this quantile of signal amplitude for gain adaptation
NORMALIZE_SIGNAL_PER_CHANNEL	1	normalize each channel (parallel signal) separately
NORMALIZE_SIGNAL_PER_CHANNEL_QUANTILE	0.1;	use this quantile of signal amplitudes for per-channel normalization
<b>ARTIFACT_TYPES</b>	{ 'ARTIF', 'UNSURE' }	cell array with artifact type abbreviations (max TODO types)
ARTIFACT_COLOR	'r'	color to plot artifact segments (shown in overview only)
ANNOT_DEFAULT_FILENAME	'sigInspectAnnot##.mat'	default name for sigInspect annotation when save button is hit. ## is replaced by current date and time
SHOW_SIG_INFO	1	show signal info from interface
overview window		
OVERVIEW_DECIMATE_FACTOR	20	decimate signals by this factor before view in the overview window - this is multiplied with the basic DECIMATE_FACTOR
OVERVIEW_GAIN	5	increase gain in overview window (to adapt for decimation)
spectrogram plot		
SPECTROGRAM_NFFT	1024	points to be used in spectrogram FFT. Window length is set to the same value, overlap is 50%
SPECTROGRAM_FREQ_LIMS	[0 3000]	limit spectrogram to this frequency range - in Hz
DISABLE_SPECTROGRAM	0	disable spectrogram (will use the free space for signal plot)
ENABLE_WHOLE_SPECTROGRAM	[0]	Adds a 'w' toggle to switch between current second and whole signal spectrogram (Beta)
audio		
START_WITH_SOUND_ON	0	start with sound turned on
USE_BEEP	0	beep when going to another signal
USE_AUDIOPAYER	1	use audio playback by audioplayer (matlab audioplayer instead of play())
PLAY_VIA_INTERNAL_PLAYER	1	use more sophisticated audio - may prevent audio dropouts and other problems on some systems
INTERNAL_PLAYER_DBG	0	show debug messages from internal player (works only if audioplayer is chosen)
PLAY_SOUND_IN_SYNC	0	synchronized audio mode
CHECK_CONST_SAMPLES	1	check signal for constant samples - visualize them in color

though the Object programming is not very common in Matlab, implementing your own interface is very simple, as you will see below.

## In OOP terms

In Object Oriented Programming (OOP) terms, such as in Java, the `sigInspectDataInterface` is an interface<sup>4</sup> you have to inherit in your own class. This is done in Matlab using the "@" sign after your function name. The `sigInspectDataInterface` has two abstract methods `getSignalIds()` and `getSignalsById(signalId)` which you have to implement to have a valid class. It may be also useful to define a constructor for your class, which allows you to e.g. pick data to view right on the object initialization. You can see an example class below.

## Reference

### getSignalIds

```
1 function signalIds = getSignalIds(obj)
```

`signalIds` is a cell array of strings - unique signal identifiers. This function is called only once during `sigInspect` initialization. After that, `sigInspect` stores the `signalIds` internally.

### getSignalsById

```
1 function [signals chInfo]= getSignalsById(obj,signalId)
```

This function is called whenever the viewed signal changes and new one has to be loaded. Here, `signals` is a  $C \times N$  matrix, where  $C$  is the number of channels and  $N$  is the number of samples. Thus, all channels have to be of the same length - use trailing zeros if one of the channels is shorter. The other output `chInfo` is a single string, which is shown in time series plot title. You can use this parameter to provide info about individual channels (such as in Figure 1, all text beyond the '18t245p3x:' (`signalId`) in the plot title.). If you do not need additional signal info, just return empty string.

### settings

The last attribute of the abstract class is the `settings` structure. By default, this is an empty structure (i.e. it is not abstract), so you do not have to define within your interface. However, it may be very useful for defining settings for your data source. Settings are defined simply by assigning the settings structure a field of the same name as the property you want to change. For example, you can change the sampling frequency within your interface as follows:

```
1 intf=myInterface();
2 intf.SAMPLING_FREQ=10000;
3 sigInspect(intf);
```

---

<sup>4</sup>in Matlab it is in fact an *abstract class* as the language defines no interfaces and ignores the subtle distinctions between the two

## Example interface

The following example code creates an interface to load signals from all csv files, present in a directory. In this example, the csv files contain three-channel micro-EEG signal in columns. The example ignores error checking (dir/file existence) for the sake of brevity.

```
1  classdef sigInspectDataCsv < sigInspectDataInterface
2      % define class, inherit sigInspectDataInterface abstract class
3      properties
4          dirPath='';
5      end
6      methods
7          % constructor - just store the path, set settings
8          function obj=sigInspectDataCsv(dirPath)
9              obj.dirPath = dirPath;
10             obj.settings.SAMPLING_FREQ=6000; % 6kHz sampling rate
11             obj.settings.PLOT_STEP=1.5;      % distance between channels on the y-axis
12             obj.settings.ARTIFACT_TYPES={'Type A','Type B','Unsure'};
13         end
14         % return list of signal ids - load all csv files from a directory,
15         % use filenames as signalIds
16         function signalIds = getSignalIds(obj)
17             lst=dir([obj.dirPath '/*.csv']);
18             signalIds = {lst(:).name}';
19         end
20
21         % read signals based on signalId (=filename)
22         function [signals chInfo]= getSignalsById(obj,signalId)
23             chInfo='';
24             signals=csvread([obj.dirPath '/' signalId]);
25         end
26     end
27 end
```

## Using the interface

Once you have your class implemented, it is very easy to use it

```
1  % initialize interface using its constructor
2  intf=sigInspectDataCsv('csvDemo/');
3  % change additional settings (optional)
4  intf.settings.NORMALIZE_SIGNAL.PER_CHANNEL = 0;
5  % run sigInspect
6  sigInspect(intf)
```

## Default interface: sigInspectDataBasic

When you provide sigInspect with a cell array of signals or a path to a mat-file with data, it will initialize the default interface: sigInspectDataBasic, which accepts these on the input (see inline help for details). The main disadvantage of calling sigInspect with signals directly is that you have no possibility to change settings, other than to hard-code them in sigInspect.m directly. A simple workaround is to initialize the basic interface yourself as follows:

```
1  % prepare data + initialize the Basic interface
2  signals = {signal1 signal2 signal3};
3  intf = sigInspectDataBasic(signals);
4  % change the settings
5  intf.settings.ARTIFACT_TYPES={'Wavy','Spiky','Other'};
6  % run sigInspect
7  sigInspect(intf);
```

## Automatic artifact annotation (Beta)

To assist you with manual artifact annotation in your signals, the sigInspect tool provides you with several methods of artifact detection. This tool is fine-tuned for microEEG data and in its initial release, it is intended as a starting-point for manual signal inspection, rather than as a stand-alone tool. The workflow is following:

1. initialize data interface of your choice
2. run sigInspectAutoLabel on this interface, which will result in automatic annotation in a mat-file.
3. run sigInspect with the same interface, load the annotation
4. update the artifact annotation manually to your taste

To see different argument types of sigInspectAutoLabel, check out its inline help. For example, your artifact labelling code could look like this:

```
1 % use the Basic interface + data in mat-file
2 intf = sigInspectDataBasic('mySignals.mat');
3 % change the settings
4 intf.settings.ARTIFACT_TYPES={'Automatic','MyArtif1','MyArtif2'};
5 % run autoLabel with the default method
6 sigInspectAutoLabel(intf,'myAutoAnnot.mat');
7 % run sigInspect
8 sigInspect(intf);
9 % now click the "load" icon in the toolbar and open myAutoAnnot.mat
```

## Artifact types handling

No matter what artifact types you set in your interface, the sigInspectAutoLabel function uses only 0-1 (clean-artifact) classification. The number of artifact types you define will only be used to resize annotation matrix for each signal to appropriate dimension. It may be a good idea to use artifact types in the form intf.ARTIFACT\_TYPES={'AUTO','myType1','myType2'}, so that you keep record of the original automatic annotation. You can also set the property intf.settings.ARTIFACT\_AUTOLABEL\_WHICH ... = 2 to set which artifact type shall be assigned by the auto label procedure. Default value is the 1.

## Artifact classification methods

While sigInspectAutoLabel is gateway function for user's convenience, artifact classification itself is performed by the function sigInspectClassify. This function also implements all classification methods mentioned below. Artifact classification is done on a set of features computed from the raw signal by sigInspectComputeFeatures - the feature set varies from method to method. Particular setting of input parameters for each method can be found in help sigInspectClassify.

### psd - deviation from ideal spectrum

The *psd* method is designed to detect artifacts with very strong frequency components and is based on our 2015 EMBC paper [2]. The classifier stores a mean Power spectrum (or Power Spectral Density - PSD), calculated on clean portions of a training set of 60 ten seconds long microEEG recordings. Manual annotation of the training data was performed by 5 trained experts using previous version of the sigInspect tool. Classification threshold was set to maximum accuracy on the training data. Detection sensitivity can be modified by adjusting detection threshold (additional input parameter): default value 0.0085, set about 0.0075 to obtain more sensitive classification (at the cost of decreased specificity). More details can be found in the paper [2] or thesis [1].



### cov - maximal stationary segment

The *cov* method is based on the solution presented by Aboy and Falkenberg in [3, 4]. The method compares two signal segments by the ratio of variance of their autocorrelation functions. SigInspect implements an extended variant presented in the paper [2], which uses the same principle to construct a similarity matrix between all possible segment pairs and returns the largest component of similarity (i.e. non-contiguous set of signal segments between which exists a sub-threshold path). Detection threshold (default: 1.2) and length of the segment (0.25s) can be set as additional parameters. Classification example using the *cov* method can be found in the following code:

```
1 % use the Basic interface + data in mat-file
2 intf = sigInspectDataBasic('mySignals.mat');
3 % change the settings
4 intf.settings.ARTIFACT.TYPES={'Automatic','Unsure'};
5 % run autoLabel with the cov method
6 sigInspectAutoLabel(intf,'myAutoAnnot.mat',[], 'cov',1.2, 0.25);
7 % sampling frequency left default (24kHz), threshold 1.2, window length 0.25s
8 % run sigInspect
9 sigInspect(intf);
10 % now click the "load" icon in the toolbar and open myAutoAnnot.mat
```

### tree - decision tree classification

The *tree* classifier is presented in the thesis [1] and uses a decision tree classifier on multiple spectral and time-domain features. The tree classifier achieved highest artifact detection accuracy on a multi-centric data set [1]. Its variant *treePrg*, trained on a single center database provides a more computationally efficient alternative, omitting features with high computational cost.

## About + license

SigInspect is distributed under the Lesser General Public License (<https://www.gnu.org/licenses/lgpl-3.0.en.html>), which means you can use sigInspect free of charge for research, educational, commercial or any other purpose, but you should keep the license when using larger portions of the code. If you publish results processed with our tool, we expect you to give us appropriate authorship credit (e.g. by citing our paper [2]).

We will be very happy if sigInspect helps you in your research. If you are happy (or unhappy) with our tool, we will be glad if you let us know!

## FAQ

**What does sigInspect do?** sigInspect is a Matlab GUI tool that Allows you to inspect and annotate multi-channel signals.

**Can I start sigInspect without signals?** No, signals or interface must be provided on startup

**Can I change interface during runtime?** No, this is currently not supported

**Signals are overlapping, what shall I do?** You can A) change signal gain using +/- keys or slider on the left or B) change the PLOT.STEP settings field

**What is the *interface*?** The term *interface* refers here to an object of class `sigInspectDataInterface`, which can be used to tailor the GUI to multiple data sources - see the "Advanced data loading" section for details.

**Do I need the *interface*?** No, you do not have to worry - just start the sigInspect the basic way (with signals as a parameter), it will initialize `sigInspectDataBasic` interface for you.

**Can I annotate segments shorter than the whole current view (1s)?** No, this is not possible. But trust us - this way it is **much** faster and it may also be really hard to identify beginning and end of an artifact.

## References

- [1] Bakstein E. *Deep Brain Recordings in Parkinson's Disease: Processing, Analysis and Fusion with Anatomical Models*, PhD thesis, Czech Technical University in Prague, 2016
- [2] Bakstein E., Schneider J., Sieger T., Novak D., Wild J. and Jech R. *Supervised Segmentation of Microelectrode Recording Artifacts Using Power Spectral Density* Proceedings of the 37th Annual International Conference of the IEEE Engineering in Medicine and Biology Society Milan, Italy, August 2015
- [3] J. H. Falkenberg and J. McNames, *Segmentation of extracellular microelectrode recordings with equal power* in Proceedings of the 25th Annual International Conference of the IEEE Engineering in Medicine and Biology Society, 2003, vol. 3, pp. 24752478.
- [4] M. Aboy and J. H. Falkenberg, *An automatic algorithm for stationary segmentation of extracellular microelectrode recordings* Med. Biol. Eng. Comput., vol. 44, no. 6, pp. 511515, Jun. 2006.